



Starting the psychopy.iohub Process ¶

To use ioHub within your PsychoPy Coder experiment script, ioHub needs to be started at the start of the experiment script. The easiest way to do this is by calling the `launchHubServer` function.

launchHubServer function

`psychopy.iohub.client.launchHubServer(**kwargs)`

The `launchHubServer` function is used to start the ioHub Process by the main psychopy experiment script.

To use ioHub for keyboard and mouse event reporting only, simply use the function in a way similar to the following:

```
from psychopy.iohub import launchHubServer

# Start the ioHub process. The return variable is what is used
# during the experiment to control the ioHub process itself,
# as well as any running ioHub devices.
io=launchHubServer()

# By default, ioHub will create Keyboard and Mouse devices and
# start monitoring for any events from these devices only.
keyboard=io.devices.keyboard
mouse=io.devices.mouse

# As a simple example, use the keyboard to have the experiment
# wait until a key is pressed.

print "Press any Key to Exit Example....."

keys = keyboard.waitForKeys()

print "Key press detected, exiting experiment."
```

`launchHubServer()` accepts several kwarg inputs, which can be used when more complex device types are being used in the experiment. Examples are eye tracker and analog input devices.

Please see the `psychopy/demos/coder/iohub/launchHub.py` demo for examples of different ways to use the `launchHubServer` function.

ioHubConnection Class

The `psychopy.iohub.ioHubConnection` object returned from the `launchHubServer` function provides methods for controlling the ioHub process and accessing ioHub devices and events.

`class psychopy.iohub.client.ioHubConnection(ioHubConfig=None, ioHubConfigAbsPath=None)`

`ioHubConnection` is responsible for creating, sending requests to, and reading replies from the ioHub Process. This class can also shut down and disconnect the ioHub Process.

The `ioHubConnection` class is also used as the interface to any ioHub Device instances that have been created so that events from the device can be monitored. These device objects can be

accessed via the `ioHubConnection .devices` attribute, providing 'dot name' attribute access, or by using the `.deviceByLabel` dictionary attribute; which stores the device names as keys,

Using the `.devices` attribute is handy if you know the name of the device to be accessed and you are sure it is actually enabled on the ioHub Process. The following is an example of accessing a device using the `.devices` attribute:

```
# get the Mouse device, named mouse
mouse=hub.devices.mouse
current_mouse_position = mouse.getPosition()

print 'current mouse position: ', current_mouse_position

# Returns something like:
# >> current mouse position: [-211.0, 371.0]
```

getDevice(deviceName)

Returns the `ioHubDeviceView` that has a matching name (based on the `device : name` property specified in the `ioHub_config.yaml` for the experiment). If no device with the given name is found, `None` is returned. Example, accessing a Keyboard device that was named 'kb'

```
keyboard = self.getDevice('kb')
kb_events= keyboard.getEvent()
```

This is the same as using the 'natural naming' approach supported by the `.devices` attribute, i.e:

```
keyboard = self.devices.kb
kb_events= keyboard.getEvent()
```

However the advantage of using `getDevice(device_name)` is that an exception is not created if you provide an invalid device name, or if the device is not enabled on the ioHub server (for example if the device hardware was not connected when the ioHub server started). Instead `None` is returned by this method. This allows for conditional checking for the existence of a requested device within the experiment script, which can be useful in some cases.

Args:

`deviceName (str)`: Name given to the ioHub Device to be returned

Returns:

`device (ioHubDeviceView)` : the PsychoPy Process representation for the device that matches the name provided.

getEvents(device_label=None, as_type='namedtuple')

Retrieve any events that have been collected by the ioHub Process from monitored devices since the last call to `getEvents()` or `clearEvents()`.

By default all events for all monitored devices are returned, with each event being represented as a `namedtuple` of all event attributes.

When events are retrieved from an event buffer, they are removed from that buffer as well.

If events are only needed from one device instead of all devices, providing a valid device name as the `device_label` argument will result in only events from that device being returned.

Events can be received in one of several object types by providing the optional `as_type` property to the method. Valid values for `as_type` are the following str values:

- ‘list’: Each event is sent from the ioHub Process as a list of ordered attributes. This is the most efficient for data transmission, but not for human readability or usability. However, if you do want events to be kept in list form, set `as_type = ‘list’`.
- ‘astuple’: Each event is converted to a namedtuple object. Event attributes are accessed using natural naming style (dot name style), or by the index of the event attribute for the event type. The namedtuple class definition is created once for each Event type at the start of the experiment, so memory overhead is almost the same as the event value list, and conversion from the event list to the namedtuple is very fast. This is the default, and normally most useful, event representation type.
- ‘dict’: Each event converted to a dict object, keys equaling the event attribute names, values being, well the attribute values for the event.
- ‘object’: Each event is converted into an instance of the ioHub DeviceEvent subclass based on the event’s type. This conversion process can take a bit of time if the number of events returned is large, and currently there is no real benefit converting events into DeviceEvent Class instances vs. the default namedtuple object type. Therefore this option should be used rarely.

Args:

`device_label (str)`: Indicates what device to retrieve events for. If `None` (the default) returns device events from all devices.

`as_type (str)`: Indicates how events should be represented when they are returned to the user. Default: ‘namedtuple’.

Returns:

tuple: A tuple of event objects, where the event object type is defined by the ‘as_type’ parameter.

clearEvents(*device_label*=‘all’)

Clears events from the ioHub Process’s Global Event Buffer (by default) so that unneeded events are not sent to the PsychoPy Process the next time `getEvents()` is called.

If `device_label` is ‘all’, (the default), then events from both the ioHub *Global Event Buffer* and all *Device Event Buffer’s* are cleared.

If `device_label` is `None` then all events in the ioHub *Global Event Buffer* are cleared, but the *Device Event Buffers* are unaffected.

If `device_label` is a str giving a valid device name, then that *Device Event Buffers* is cleared, but the *Global Event Buffer* is not affected.

Args:

`device_label (str)`: device name, ‘all’, or `None`

Returns:

`None`

sendMessageEvent(*text*, *category*=“, *offset*=0.0, *sec_time*=None)

Create and send an Experiment MessageEvent to the ioHub Server Process for storage with the rest of the event data being recorded in the ioDataStore.

Note

MessageEvents can be thought of as DeviceEvents from the virtual PsychoPy Process “Device”.

Args:

text (str): The text message for the message event. Can be up to 128 characters in length.

category (str): A 0 – 32 character grouping code for the message that can be used to sort or group messages by ‘types’ during post hoc analysis.

offset (float): The sec.msec offset to apply to the time stamp of the message event. If you send the event before or after the time the event actually occurred, and you know what the offset value is, you can provide it here and it will be applied to the ioHub time stamp for the MessageEvent.

sec_time (float): The time stamp to use for the message in sec.msec format. If not provided, or None, then the MessageEvent is time stamped when this method is called using the global timer.

Returns:

bool: True

createTrialHandlerRecordTable(trials)

Create a condition variable table in the ioHub data file based on the a psychopy TrialHandler. By doing so, the ioHub data file can contain the DV and IV values used for each trial of an experiment session, along with all the ioHub device events recorded by ioHub during the session. Example psychopy code usage:

```
# Load a trial handler and
# create an associated table in the ioHub data file
#
from psychopy.data import TrialHandler, importConditions

exp_conditions=importConditions('trial_conditions.xlsx')
trials = TrialHandler(exp_conditions,1)

# Inform the ioHub server about the TrialHandler
#
io.createTrialHandlerRecordTable(trials)

# Read a row of the trial handler for
# each trial of your experiment
#
for trial in trials:
    # do whatever...

# During the trial, trial variable values can be updated
#
trial['TRIAL_START']=flip_time

# At the end of each trial, before getting
# the next trial handler row, send the trial
# variable states to ioHub so they can be stored for future
# reference.
#
io.addTrialHandlerRecord(trial.values())
```

addTrialHandlerRecord(cv_row)

Adds the values from a TrialHandler row / record to the ioHub data file for future data analysis use.

Parameters: cv_row –

Returns: None

getTime()

Deprecated Method: Use `Computer.getTime` instead. Remains here for testing time bases between processes only.

setPriority(*level*='normal', *disable_gc*=False)

See `Computer.setPriority` documentation, where current process will be the iohub process.

getPriority()

See `Computer.getPriority` documentation, where current process will be the iohub process.

enableHighPriority(*disable_gc*=False)

Deprecated Method: Use `setPriority('high', disable_gc)` instead.

disableHighPriority()

Deprecated Method: Use `setPriority('normal')` instead.

enableRealTimePriority(*disable_gc*=False)

Deprecated Method: Use `setPriority('realtime', disable_gc)` instead.

disableRealTimePriority()

Deprecated Method: Use `setPriority('normal')` instead.

getProcessAffinity()

Returns the current **ioHub Process Affinity** setting, as a list of 'processor' id's (from 0 to `getSystemProcessorCount()-1`). A Process's Affinity determines which CPU's or CPU cores a process can run on. By default the ioHub Process can run on any CPU or CPU core.

This method is not supported on OS X at this time.

Args:

None

Returns:

list: A list of integer values between 0 and `Computer.getSystemProcessorCount()-1`, where values in the list indicate processing unit indexes that the ioHub process is able to run on.

setProcessAffinity(*processor_list*)

Sets the **ioHub Process Affinity** based on the value of `processor_list`. A Process's Affinity determines which CPU's or CPU cores a process can run on. By default the ioHub Process can run on any CPU or CPU core.

The `processor_list` argument must be a list of 'processor' id's; integers in the range of 0 to `Computer.processing_unit_count-1`, representing the processing unit indexes that the ioHub Server should be allowed to run on. If `processor_list` is given as an empty list, the ioHub Process will be able to run on any processing unit on the computer.

This method is not supported on OS X at this time.

Args:

`processor_list` (list): A list of integer values between 0 and `Computer.processing_unit_count-1`, where values in the list indicate processing unit indexes that the ioHub process is able to run on.

Returns:

None

flushDataStoreFile()

Manually tell the ioDataStore to flush any events it has buffered in memory to disk.”

Args:

None

Returns:

None

shutdown()

Tells the ioHub Process to close all ioHub Devices, the ioDataStore, and the connection monitor between the PsychoPy and ioHub Processes. Then exit the Server Process itself.

Args:

None

Returns:

None

quit()

Same as the shutdown() method, but has same name as PsychoPy core.quit() so maybe easier to remember.